


More flexibility*: CM Use for CT motion™ and Max 3™

Make working easier by using different contrast media for different patients. Flexible contrast media management allows you to optimize your workflow, produce less waste and work efficiently.*




Tandem / CM-Select function

- Changeover of CM with the same active ingredients at the touch of a button
- The 24h pump hose does not need to be changed / an Easy-Click-Cassette is required
- Automatic rinsing program

Extended Tandem / CM-Select function

- Use of CM with different active ingredients
- Changeover of contrast media using a special rinsing program
- The 24h pump hose does not need to be changed / an Easy-Click-Cassette is required



Seamless / CM-Loop function

- Double CM reserve thanks to automatic switchover with identical CM
- Low quantities of discarded CM
- Efficiency thanks to high patient throughput and patient change-overs in a short space of time

Let us help you boost your efficiency:

- Activation of the Seamless and Tandem functions
- User training tailored to the functions and in line with your requirements
- Please do not hesitate to contact your sales representative for more detailed information

* Compared to our devices without (Extended) Tandem / CM-Select function.

Contact us:

