

uCentum™ components

Crosslinks CS 3810 - [15...20] -S / 6 products

- Clamping of rods with Ø 5.5 mm and Ø 6.0 mm
- Six different sizes (XXS, XS, S, M, L, XL)

Overview



CS 3810-15-S

Crosslink variable, for Ø 5.5 mm/Ø 6.0 mm, sterile, XXS, length 35-43 mm



CS 3810-16-S

Crosslink variable, for Ø 5.5 mm/Ø 6.0 mm, sterile, XS, length 41-56 mm



CS 3810-17-S

Crosslink variable, for Ø 5.5 mm/Ø 6.0 mm, sterile, S, length 53-68 mm



CS 3810-18-S

Crosslink variable, for Ø 5.5 mm/Ø 6.0 mm, sterile, M, length 66-80 mm



CS 3810-19-S


Crosslink variable, for Ø 5.5 mm/Ø 6.0 mm, sterile, L, length 78-94 mm



CS 3810-20-S

Crosslink variable, for Ø 5.5 mm/Ø 6.0 mm, sterile, XL, length 92-107 mm

Determining the correct implant size



The distance between the two rods can be determined using the measuring caliper (CS 5788), **taking the measurement on the outside of the rods.**

dedicated to you.

Ulrich
medical