


uCentumTM

comprehensive posterior system

uBaseTM

iliosacral fixation system




Complete
system for open
and percutaneous
applications

uCentum™

System

- Posterior stabilization of the thoracic to lumbosacral spine
- Complete system for open and percutaneous applications
- Ø 5.5 mm and Ø. 6.0 mm rods
- Polyaxial and monoaxial screws
- Pre-fixation screws and long-arm screws
- MIS X-Tab screw for percutaneous application (also pre-fixation)
- Can be combined with all the screws from the uBase system
- Use in combination with the Brainlab navigation software possible (both open and percutaneous)
- Extensive navigation module for preparation of the surgical procedure and screw insertion
- Combination with the cervicothoracic screw-rod system neon³
- Implants available sterile



Cannulated and perforated screws (except Ø 4.5 mm)




MIS X-Tab screw for percutaneous application



Offset connector

Low-profile screw head

High flexibility due to various
connector lengths and additional
connection options



uBase™

System

- Stabilization of the sacrum and ilium in connection with uCentum
- Purely sacral, sacroiliac or purely iliac placement
- Polyaxial screw head, pre-angled medial/lateral or cranial/caudal
- Additionally available: fixed-angle screw
- All screws cannulated
- Implants available sterile

Advantages


- Improved safety thanks to cannulated screws and extensive navigation module for preparation of the surgical procedure and screw insertion
- 48° polyaxiality
- Percutaneous application (polyaxial screws)
- Screw head with low profile




Iliosacral screw, polyaxial



Purely sacral screw placement



Sacroiliac screw placement



Purely iliac screw placement



uCentum™

Advantages

- Improved safety thanks to cannulated screws and extensive navigation module for preparation of the surgical procedure and screw insertion
- Cement augmentation through the screw thanks to lateral perforations
- Rods with and crosslinks for Ø 5.5 and Ø 6.0 mm
- Combination with the cervicothoracic screw-rod system neon³
- New reduction instrumentation



uCentum™

Components

Implants

Screw, polyaxial, Ø 4.5 mm, length 25 mm – 50 mm
Screw, polyaxial, Ø 5.5 mm, length 25 mm – 55 mm, perforated
Screw, polyaxial, Ø 6.5 mm, length 30 mm – 60 mm, perforated
Screw, polyaxial, Ø 7.5 mm, length 35 mm – 60 mm, perforated
Screw, polyaxial, Ø 8.5 mm, length 35 mm – 60 mm, perforated
Screw, polyaxial, Ø 10 mm, length 45 mm – 60 mm, perforated
Screw, polyaxial, sterile, Ø 11.5 mm, length 45 mm – 55 mm, perforated*

Product no.

CS 3802-045-025 – CS 3802-045-050
CS 3802-055-025 – CS 3802-055-055
CS 3802-065-030 – CS 3802-065-060
CS 3802-075-035 – CS 3802-075-060
CS 3802-085-035 – CS 3802-085-060
CS 3802-100-045 – CS 3802-100-060
CS 3802-115-045-S – CS 3802-115-055-S

Screw, monoaxial, Ø 4.5 mm, length 25 mm – 50 mm
Screw, monoaxial, Ø 5.5 mm, length 25 mm – 55 mm, perforated
Screw, monoaxial, Ø 6.5 mm, length 30 mm – 60 mm, perforated
Screw, monoaxial, Ø 7.5 mm, length 35 mm – 60 mm, perforated
Screw, monoaxial, Ø 8.5 mm, length 35 mm – 60 mm, perforated
Screw, monoaxial, Ø 10 mm, length 45 mm – 60 mm, perforated

CS 3803-045-025 – CS 3803-045-050
CS 3803-055-025 – CS 3803-055-055
CS 3803-065-030 – CS 3803-065-060
CS 3803-075-035 – CS 3803-075-060
CS 3803-085-035 – CS 3803-085-060
CS 3803-100-045 – CS 3803-100-060

Screw, pre-fixation, Ø 4.5 mm, length 25 mm – 50 mm
Screw, pre-fixation, Ø 5.5 mm, length 25 mm – 55 mm, perforated
Screw, pre-fixation, Ø 6.5 mm, length 30 mm – 60 mm, perforated
Screw, pre-fixation, Ø 7.5 mm, length 35 mm – 60 mm, perforated
Screw, pre-fixation, Ø 8.5 mm, length 35 mm – 60 mm, perforated
Screw, pre-fixation, Ø 10 mm, length 45 mm – 60 mm, perforated

CS 3804-045-025 – CS 3804-045-050
CS 3804-055-025 – CS 3804-055-055
CS 3804-065-030 – CS 3804-065-060
CS 3804-075-035 – CS 3804-075-060
CS 3804-085-035 – CS 3804-085-060
CS 3804-100-045 – CS 3804-100-060

Long-arm screw, pre-fixation, Ø 5.5 mm, length 25 mm – 55 mm, perforated
Long-arm screw, pre-fixation, Ø 6.5 mm, length 30 mm – 60 mm, perforated
Long-arm screw, pre-fixation, Ø 7.5 mm, length 35 mm – 60 mm, perforated
Long-arm screw, pre-fixation, Ø 8.5 mm, length 35 mm – 60 mm, perforated

CS 3805-055-025 – CS 3805-055-055
CS 3805-065-030 – CS 3805-065-060
CS 3805-075-035 – CS 3805-075-060
CS 3805-085-035 – CS 3805-085-060

Screw, MIS X-Tab, polyaxial, sterile, Ø 4.5 mm, length 25 mm – 50 mm*
Screw, MIS X-Tab, pre-fixation, sterile, Ø 4.5 mm, length 25 mm – 50 mm*
Screw, MIS X-Tab, polyaxial, sterile, Ø 5.5 mm, length 25 mm – 55 mm, perforated*
Screw, MIS X-Tab, pre-fixation, sterile, Ø 5.5 mm, length 25 mm – 55 mm, perforated*
Screw, MIS X-Tab, polyaxial, sterile, Ø 6.5 mm, length 30 mm – 55 mm, perforated*
Screw, MIS X-Tab, pre-fixation, sterile, Ø 6.5 mm, length 30 mm – 55 mm, perforated*
Screw, MIS X-Tab, polyaxial, sterile, Ø 7.5 mm, length 35 mm – 55 mm, perforated*
Screw, MIS X-Tab, pre-fixation, sterile, Ø 7.5 mm, length 35 mm – 55 mm, perforated*

CS 3806-045-025-S – CS 3806-045-050-S
CS 3806-045-925-S – CS 3806-045-950-S
CS 3806-055-025-S – CS 3806-055-055-S
CS 3806-055-925-S – CS 3806-055-955-S
CS 3806-065-030-S – CS 3806-065-055-S
CS 3806-065-930-S – CS 3806-065-955-S
CS 3806-075-035-S – CS 3806-075-055-S
CS 3806-075-935-S – CS 3806-075-955-S

Combination with uBase

Iliosacral screw, Ø 7.5 mm, length 40 mm – 80 mm
Iliosacral screw, Ø 8.5 mm, length 40 mm – 120 mm
Iliosacral screw, Ø 10.0 mm, length 80 mm – 140 mm

CS 7170-075-040 – CS 7170-075-080
CS 7170-085-040 – CS 7170-085-120
CS 7170-100-080 – CS 7170-100-140

Iliosacral screw, polyaxial, Ø 7.5 mm, length 40 mm – 80 mm, sterile*
Iliosacral screw, polyaxial, Ø 8.5 mm, length 40 mm – 120 mm, sterile*
Iliosacral screw, polyaxial, Ø 10.0 mm, length 80 mm – 140 mm, sterile*

CS 7172-075-040-S – CS 7172-075-080-S
CS 7172-085-040-S – CS 7172-085-120-S
CS 7172-100-080-S – CS 7172-100-140-S

Iliac screw, polyaxial, sterile, Ø 8.5 mm, length 80 mm – 120 mm
Iliac screw, polyaxial, sterile, Ø 10.0 mm, length 80 mm – 120 mm

CS 7173-085-080-S – CS 7173-085-120-S
CS 7173-100-080-S – CS 7173-100-120-S

*These products are only available sterile

Information on sterile packed implants

To order sterile packed implants, add "S" to the end of the product number (example: CS 3802-045-025-S).



ulrich GmbH & Co. KG
Buchbrunnenweg 12
89081 Ulm, Germany
T +49 731 9654-0
spine@ulrichmedical.com
www.ulrichmedical.com

ulrich medical France SAS
10, Allée Aristide Maillol
31770 Colomiers, France
T +33 5 34 50 91 02
info@ulrichmedical.fr
www.ulrichmedical.fr